

Taking a test



WAKE UP: The device will wake up and alert you with several beeps when you are required to take a test.



BLOW: The device is ready to accept a breath test.
PICK UP THE ENTIRE Unit – HANDSET AND BASE.
 Blow into the device for about 5 seconds until you hear a beep. The screen will say "WAIT".



PASS: You have passed the test and the screen will display your next test time. This could be an exact time or a time range.



FAILED: You have failed the test. The screen will say "WAIT" while the device resets to allow you to take a retest. The screen will count down the number of seconds you have to wait.





NOTE: Trace amounts of alcohol can give a positive reading. To avoid this, follow these guidelines:

- Don't eat or drink five minutes before testing
- Rinse your mouth with water before taking a test
- Ensure proper ventilation

May Contain Alcohol (this list is not all inclusive)

- Hand Sanitizer
- Antibacterial Soap
- Mouthwash
- Hot Sauce
- Sugarless Gum
- Household Products (air fresheners, bleach, laundry detergent, glue)
- Hygiene Products (deodorant sprays, aftershave, perfume, toothpaste, insect repellent)
- Medicines and Inhalers
- Energy Drinks
- Fruit, Bread, Pizza
- E-cigarettes

RETEST: The screen will display **BLOW**. The device is ready to accept a breath test*. The screen will display **FAILED** or **PASS** depending on the outcome of the test. The Screen will display your next test time.



*The number of and /or time allowed for retests will vary by jurisdiction and is determined by your monitoring authority.

Required Monitoring Appointments

When you receive your portable alcohol monitoring device, a service appointment will be made for you. This is typically 30 days or less from the date of receipt and every 30 days thereafter until the end of your program.

During this appointment the device's memory will be downloaded and analyzed. Results of breath tests, violations and other data will be reported to your monitoring authority. The service provider will inspect the device for signs of any attempted tampering and, if discovered, this will be reported as well. The device will also be calibrated, inspected for proper functioning and your next payment will be collected.

You must call **24 hours in advance** to change a scheduled monitoring appointment. Failure to do so will result in a Missed Appointment Fee.

WARNING: If you miss a scheduled monitoring appointment or fail to respond to an early recall, a non-compliance report will be sent to authorities.



Quick Reference Guide

This guide explains the use of the LifeSafer Portable Alcohol Monitoring device (PAM). You will be trained on proper use of the LifeSafer PAM device. A schedule of test times will be established by your monitoring authority and provided with this handbook. You may also be required to take random tests throughout the day. You will be required to have the device serviced periodically. A date will be programmed into the device telling you when you are scheduled to return for service.

How it Works

The LifeSafer PAM is a breath analyzer linked to a camera. The handset features a display screen and keypad for easy navigation.



24-Hour Service Assistance

If you have any questions regarding your device, let us know! We are available to help you 24 hours a day, 365 days a year.

Our helpful customer care agents can also assist in the following:

- Scheduling appointments
- Making payments
- Troubleshooting
- Device service notifications

If you call after business hours, please listen carefully to the prompts which will allow you to leave a message, or speak to a customer care agent.

LifeSafer
Portable Alcohol Monitoring

For more information about LifeSafer, visit us online at lifesaffer.com or call us

(855) 527-1597



Portable Alcohol Monitoring Unit

Quick Reference Guide



LifeSafer
Portable Alcohol Monitoring

lifesaffer.com